

[Home](#) > [Digital Imaging](#) > [RVG Imaging](#) > RVG 6500

RVG 6500 System*



The RVG 6500 System features the latest innovations in digital radiography—while still delivering the highest image resolution (> 20 lp/mm). Since it's Wi-Fi enabled, you can also keep your operatories cable-free. The result: faster, better mobile performance with no loss in image quality. With a system this flexible and affordable, the RVG 6500 System is the ideal solution for even the most demanding applications.

Special promotion for the RVG 6500 System:**
call 1.800.944.6365 for details.

Contact Us

Digital Imaging & Software:

- 1.800.944.6365 (US)
- 1.855.762.2291 (Canada)

[Email Us](#)

Film and Anesthetics:

- 1.800.933.8031

[Email Us](#)
[Request a Demo](#)
[Email Signup](#)
[Features and Benefits](#)
[Image Gallery](#)
[Resource Library](#)
[Testimonials](#)

Resource Library

Click on the desired image or text below to view the selected subject area

[Documentation \(US\)](#)
[Documentation \(Canada\)](#)
[Technical Specifications](#)

Technical Specifications

Over the years, our name has become synonymous with providing high-quality imaging technology to the dental industry. With our RVG Systems, we offer three solutions to help meet the unique needs of your practice.

RVG 6500 System's Sensor Features	Size 0	Size 1	Size 2
True Image Resolution	15 lp/mm	Greater than 20 lp/mm	Greater than 20 lp/mm
Pixel Size	18.5 µm	18.5 µm	18.5 µm
Outside Dimensions	22.2 x 30.8 mm	27.5 x 37.7 mm	32.2 x 44.1 mm
Dimensions of Active Area	17 x 22 mm	22 x 30 mm	27 x 36 mm
Matrix Dimensions (Pixels)	900 x 1200 pixels	1200 x 1600 pixels	1440 x 1920 pixels
Battery	Lithium		
Battery Life	Approx. 180 shots		
Control Box Dimensions	Non-IPS: 83 (H) x 47 (W) x 16 (D) mm With IPS: 83 (H) x 47 (W) x 21 (D) mm		
Wireless Technology	Wi-Fi 8012-11g		

*Not available in all countries. Please check with your sales representative.

**Please note that all special promotions apply to U.S. customers only.

This website makes use of cookies to help us improve the website and give you the best experience. By using the website you agree to our use of cookies. [Read More](#)