


S-Max SG20 Contra-angle 20:1 Reduction

- Double Sealing System
- Stressless Body Shape
- Excellent Durability

SG50MS Micromotor

- Compact, Short and Lightweight Micromotor
- High Power 210W / Powerful Torque 50Ncm
- Excellent Durability and Autoclavability



Double Sealing System



Surgic AP Complete Set Non-Optic

Surgic AP Non-Optic Complete Set With S-Max SG20 Handpiece

MODEL: **Surgic AP (230V)** ORDER CODE: **Y1002184**

Contents:

- Control Unit • Non-Optic Micromotor (SG50MS)
- SG20 Handpiece (20:1 Reduction)
- FC-55 Foot Control (PXB certified) • Handpiece Stand
- Irrigation Tube (5 pcs.) and Other Accessories



Specifications

Control Unit

- Max. Output/Max. Torque : 210W / 50 Ncm
- Power Supply : AC230V 50/60 Hz
- Max. Pump Output : 75 mL/min
- Programs : 8 Programs
- Dimensions : W276 x D207 x H101 (mm)

Micromotor (SG50MS)

- Motor Speed : 200-40,000 min⁻¹
- Non-optic
- Brushless

Foot Control (FC-55)

- Speed Control: Variable
- PXB certified

Carrying Case (Optional)

The NSK Carrying Case can accommodate all Surgic AP components as well as the optional sterilization cassette.



MODEL Carrying Case (Surgic AP)

ORDER CODE **Y1002197**

- Dimensions: W534 x D427 x H207 (mm)

Spare Parts

MODEL	PRODUCT	ORDER CODE	MODEL	PRODUCT	ORDER CODE	MODEL	PRODUCT	ORDER CODE
SG50MS		E290009	Irrigation Pump (Pack of 2)		Y1002196	Tube Holder (Pack of 7)		Y900083
Irrigation Tube (Pack of 5)		Y900197	Protection Plug		10001595			